



Chlorine tablets for powerful disinfection.
Effective against bacteria, fungi, viruses and spores.

terralin® NaDCC tablets

Our Plus

- Effervescent chlorine releasing tablets
- Broad spectrum of efficacy including bacteria, viruses, fungi and spores
- Convenient and easy to use, just add water
- economical
- Easy to store and distribute, no spillages or leaks

Instructions for Dilution

Disinfection of	Required Conc. Available Chlorine	Dilution Rate	Instructions
Blood spillage	10000 (P.P.M.)	18 tablets in 500 ml of water	Pour over blood, using gloves, wipe up with disinfectant saturated disposable cloth
Work surfaces, cupboards, floors, etc.	140 (P.P.M.)	1 tablet in 2 Ltr. of water	Wash down
Soiled linen	140 (P.P.M.)	1 tablet in 2 Ltr. of water	Immerse for 1 hour prior to washing
Utensils porcelain, glassware, rubber and plastic tubing	140 (P.P.M.)	1 tablet in 2 Ltr. of water	Immerse for 1 hour before washing
General use	1000 (P.P.M.)	1 tablet in 280ml of water	In compliance with Howie Code
Pipette jars	2500 (P.P.M.)	9 tablet in 1 Ltr. of water	Drop tablets into water-filled pipette jar, Discard daily
Dishcloths, mops, etc.	60 (P.P.M.)	1 tablet in 4.6Ltr. of water	Soak to bleach clean and deodorise

Do not immerse animal fibres such as wool or silk due to strong proteolytic action.
Avoid prolonged contact with stainless steel items.

Application areas

Available Chlorine 300mg

Solutions prepared from terralin™ NaDCC tablets are suitable for the disinfection of working surfaces, blood spillage, floors, walls, utensils, linen, glassware and equipments. For specific dilution and usage, refer to "Instructions For Dilution" table.

Information for order

Item	Delivery form	Item no.
terralin NaDCC tablets	20/Carton	70002239

Product data

Composition: Sodium Dichloroisocyanurate 500mg. Excipients Q.S

Special advice

Use biocides safely. Always read the label and product information before use. For all of the related safety information please refer to the Safety Data Sheet.

Wear Protective clothing and gloves. Solution may bleach fabric. For professional use only. Contact with acids liberates toxic gas. Wash hands thoroughly after handling. Store in a well ventilated place. Keep container tightly closed. Dispose of contents/container in accordance with local regulations.

Environmental information

Schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality.



Schülke & Mayr GmbH is certified according to DIN EN ISO 9001, DIN EN ISO 14001 and DIN EN ISO 13485 (Reg.-No. 004567-MP23) and has a validated environmental system in accordance with the Eco Audit Regulation (Reg.-No. DE-150-00003).

Schülke & Mayr GmbH
22840 Norderstedt, Germany
Phone +49 (0) 40 - 521 00 - 0
www.schuelke.com
mail@schuelke.com

Schulke India Pvt. Ltd.
A-24/10, Mohan Co-op. Ind. Est.,
Mathura Road, New Delhi-110044.
Phone +91 - 11 - 4978 6666
customercare.india@schuelke.com